

A0028CC Silver Line

Typical analysis

Phosphate-free, alkaline

Recommended for cleaning contaminated

laboratory glass and plasticware

Recommended dilution 2% to 20%

Packaging

5l clear plastic bottle

Physical properties, composition and data

Mixture

Alkaline base containing anionic and non-ionic surfactants

Weight per litre

1.06kg

Laboratory preparation, applications and practices

Usage

Analytical reagent. Solvent. Extraction processes

Filter paper

Filtech No. 8613, 8614

Storage and handling information

Storage

Store in a cool place

Safety phrases

26

Risk phrases

36/38

Disposal methods

30

Transport regulations

Tariff code

3402.20.90

Hazardous class

8

Packing group

II

UN number

1824

