

A0037NN Gold Line

egg powder

Typical Analysis

Assay of albumin (ex N)	80%
Ash	3%
Loss on drying (105°C)	5%
Water (H ₂ O)	8%

Store in refrigerator

Packaging

100g white plastic jar

250g white plastic jar

Physical properties, composition and data

Solid mixture	Contains 75% ovalbumin, ovoconalbumin, ovomucid, ovomucin, ovoglobulin, B lysozyme and avadin
Synonym	Ovalbumin
Appearance	Creamy, cloggy powder
Odour	Musty odour
Solubility	Dissolves slowly in water
pH (aqueous solution)	6.8 - 7.2
Products of decomposition	Decomposes in moist air

Laboratory preparation, applications and practices

Usage	Analytical reagent. Activating agent for certain enzymes
Filter paper	Filtech no: 0222, 0225, 1839

Storage and handling information

Storage	Store in a cool place, protected from direct sunlight and moisture
Disposal methods	16

Transport regulations

Tariff code	3502.10.91
-------------	------------

