


A0039NN

Gold Line

**Packaging**

500ml amber glass bottle

25 litre plastic drum

**Physical properties, composition and data**

Chemical formula	$C_{25}H_{54}ClN$
Synonym	Methyltrioctylammonium chloride; Tricaprylmethylammonium chloride
Mixture	A fatty quaternary ammonium chloride being methyl tricaprylyl ammonium chloride stable in both acidic and alkaline media
Atomic weight	142.24
Appearance	Colourless liquid
Odour	Characteristic odour
Miscibility	Miscible with most organic solvents
Melting point	-46°C
Boiling point	168.1°C
Flash point	49°C
Ignition point	345°C

**Laboratory preparation, applications and practices**

Usage	Analytical reagent. phase transfer catalyst. Metal extraction reagent
Filter paper	Filtech no: 0165, 0225, 1839, 1845

**Storage and handling information**

Storage	Store in a cool place
Disposal methods	16

**Transport regulations**

Tariff code	3822.00.00
-------------	------------