


**A0041NN Gold Line****Typical Analysis**

Colour index 42750

**Packaging**

25g black plastic jar

**Physical properties, composition and data**Chemical formula  $C_{32}H_{28}N_3NaO_4S$ 

573.65

Atomic weight

Appearance

Dark blue powder

pH (aqueous solution)

11.5 - 13.0

**Laboratory preparation, applications and practices**

Usage

Microscopy. pH indicator. Printing inks

**Storage and handling information**

Storage

Keep container tightly closed

Disposal methods

16

**Transport regulations**

Tariff code

3204.12.00

Caution

Possible mutagen

