

A0140NN Gold Line

powder

Typical analysis

Colour index 42775

Packaging

25g black plastic jar

A0141NN Gold Line


indicator solution

Typical analysis

Colour index 42775

Packaging

100ml amber glass bottle



Physical properties, composition and data

Chemical formula	$C_{32}H_{26}N_3Na_2O_9S$
Chemical group	Stains
Synonym	Acid blue 22; Porrier's blue; Water blue
Atomic weight	737.74
Appearance	Brownish-purple to violet powder; blue liquid
Solubility/miscibility	Soluble in water and alcohol

Laboratory preparation, applications and practices

Usage Analytical reagent. Microscopy. Collagen staining

Storage and handling information

Disposal methods 16

Transport regulations

Tariff code 3219.90.90

