

# bicarbonate indicator solution

**B0238NN** **Gold Line**

## Typical analysis

Colour change colourless to red  
Sensitivity to CO<sub>2</sub> passes test

## Packaging

500ml clear plastic bottle

## Physical properties, composition and data

Chemical group Salts  
Mixture Phenolphthalein and sodium carbonate  
Appearance Dark pinkish-red liquid

## Laboratory preparation, applications and practices

Laboratory preparation Mix 20ml 1% phenolphthalein solution with 10ml 0.05M sodium carbonate solution and dilute to 100ml with distilled water  
Usage Titrations  
Filter paper Filtech no: 0222, 0225, 1839

## Transport regulations

Tariff code 3822.00.00

