

B0252NN Gold Line**Typical analysis**

Colour Index 51010

Packaging

25g black plastic jar

Physical properties, composition and data

Chemical formula	$C_{17}H_{20}ClN_3O$
Chemical group	Stains
Synonym	1,3-Diamino-7-(diethylamino)-8-methylphenoxazin-5-ium chloride
Mixture	Oxazine mixture of aminodimethyl aminomethyl diphenoxazonium chlorides
Atomic weight	317.79
Appearance	Green crystalline powder
Solubility	Soluble in water
Melting point	195 - 197°C

Laboratory preparation, applications and practices

Usage Microscopy staining

Storage and handling informationRisk phrases 20/21/22
Disposal methods 16**Transport regulations**Tariff code 3204.13.00
ECB number 233-358-3