

C0534NN Gold Line

pentadecahydrate

Typical analysis

Assay	98%
Chloride (Cl)	0.02%
Chromate (CrO ₄)	0.05%
Iron (Fe)	0.02%



Packaging

- 100g white plastic jar
- 500g white plastic jar

Physical properties, composition, and data

Chemical formula	Cr ₂ (SO ₄) ₃ · 15H ₂ O or Cr ₂ O ₁₂ S ₃ · 15H ₂ O
Synonym	Chromium III sulphate
Atomic weight	622.41
Appearance	Dark green granular powder or flakes
Solubility	Soluble in water and alcohol
Density (g/ml)	1.5 - 2.5

Laboratory preparation, applications and practices

Usage	Analytical reagent
Filter paper	Filtech no: no: 0285, 1803

Storage and handling information

Disposal methods	3
------------------	---

Transport regulations

Tariff number	2833.23.00
---------------	------------

