

chromium

C0536NN Platinum Line

1000ppm solution

Warranty certificate

Solution containing 1000ppm chromium w/v in approximately 1N HNO₃. Metal content guaranteed within 0.5% of nominal value

Packaging

500ml clear plastic bottle

Physical properties, composition and data

Mixture	Nitric acid and chromium
Weight per litre	1.023kg
Appearance	Orange liquid

Laboratory preparation, applications and practices

Usage	Analytical reagent
Filter paper	Filtech no: no: 0222, 0225, 1839

Storage and handling information

Safety phrases	26
Risk phrases	36/38
Disposal methods	1

Transport regulations

Tariff number	3822.00.00
UN number	1760

