

D0728NN Gold Line

DPD **comparator** tablets PALIN No. 1

Packaging

100 tablets

D0729NN Gold Line

DPD **comparator** tablets PALIN No 2

Packaging

100 tablets

D0730NN Gold Line

DPD **comparator** tablets PALIN No 3

Packaging

100 tablets

D0731NN Gold Line

DPD **comparator** tablets PALIN No 4

Packaging

100 tablets

D0732NN Gold Line

DPD **nessleriser** tablets PALIN No. 1

Packaging

100 tablets

D0733NN Gold Line

DPD **nessleriser** tablets PALIN No. 2

Packaging

100 tablets

D0734NN Gold Line

DPD **nessleriser** tablets PALIN No 3

Packaging


100 tablets

D0735NN Gold Line

DPD **nessleriser** tablets PALIN No 4

Packaging

100 tablets



Usage

Determination of free chlorine, or in conjunction with No. 3, determines total residual chlorine.

Usage

Determination of mono chloramine by subtracting the value (reading) of No. 1

Usage

Determination of dichloramine by subtracting the value (reading) of No. 2, or in conjunction with No. 1, determines total residual chlorine.

Usage

Determination of total residual chlorine.