

# dpx mountant

**D0738NN**

**Gold Line**

## Packaging

100ml amber glass bottle

500ml amber glass bottle

## Physical properties, composition, and data

Liquid mixture/solution	Polystyrene solution in xylene
Appearance	Water white
Melting point	80 - 88°C
Flash point	24°C
Ignition point	465°C
Density (g/ml)	0.94 at 25°C
Hazardous material	Flammable
Other information	Not recommended for thionin or iron stains

## Laboratory preparation, applications and practices

Laboratory preparation	Romanowsky, haematoxylin, eosin-methyl blue, Masson and most preparations where Canada balsam is normally used
Usage	Mounting media for microscopy

## Storage and handling information

Safety phrases	61-10-20/21-38-62
Risk phrases	53-36/37-45



170