


E0740NN	Platinum Line
powder	
Warranty certificate	
Colour Index	45380
Dye content	min. 70%
Loss on drying	10%

Packaging
25g black plastic jar
100g black plastic jar

E0741NN	Gold Line
indicator solution	
Typical analysis	
Colour Index	45380

Packaging
100ml amber glass bottle



Physical properties, composition and data

Chemical formula	Powder: $C_{20}H_6Br_4Na_2O_5$ and Solution: $C_{20}H_{14}O_5Br_4Na_2$
Chemical group	Stains
Synonym	Bromoesin; Bromofluoresceic acid; Eosin; Eosin Y; Tetrabromofluorescein
Atomic weight	691.86
Weight per litre	Solution: 1.02kg
Appearance	Red crystals with bluish tinge or brownish red powder. Red aqueous solution
Solubility/miscibility	Water
pH (aqueous solution)	9.2
Density (g/ml)	0.71

Laboratory preparation, applications and practices

Laboratory preparation	Bromination of fluorescein
Usage	Indicator
Filter paper	Filtech no: 0222, 0225, 1839

Storage and handling information

Safety phrases	22-24
Risk phrases	20/22
Disposal methods	16

Transport regulations

Tariff code	3204.12.00
-------------	------------

