


E0745NN**Platinum Line****Warranty certificate**

Colour Index	15705
Loss on drying	10%
Sensitivity as complexometric indicator	passes test

Packaging

25g black plastic jar

**Physical properties, composition and data**

Chemical formula	$C_{20}H_{13}N_3NaO_5S$
Chemical group	Azo dyes
Atomic weight	416.39
Appearance	Black powder
Solubility	Soluble in water making a purplish colour

Laboratory preparation, applications and practices

Usage	Complexometric indicator
-------	--------------------------

Storage and handling information

Risk phrases	36-37-38
Disposal methods	16

Transport regulations

Tariff code	3204.19.00
-------------	------------

