

F0870NN Gold LineTypical analysis
Colour Index 42053**Packaging**

5g black plastic jar

Physical properties, composition and data

Chemical formula	$C_{37}H_{34}N_2Na_2O_{10}S_3$
Chemical group	Stains
Synonym	Food green 3
Atomic weight	808.85
Appearance	Red brown powder
Solubility	Soluble in water, ethanol, dimethylsulphoxide
Melting point	Decomposes
pH (aqueous solution)	2.7

Laboratory preparation, applications and practices

Usage Biological stain

Storage and handling information

Risk phrases	36/37/38
Disposal methods	16

Transport regulations

Tariff code 3822.00.00

