

fehling's solution "a"

F0875NN **Gold Line**

Typical analysis

Assay 69.12 - 69.40%

Packaging

500ml clear plastic bottle

Physical properties, composition and data

Liquid mixture/solution	A copper sulphate solution
Density (g/ml)	1.037 - 1.043
Appearance	Clear blue liquid
Miscibility	Miscible with water

Laboratory preparation, applications and practices

Laboratory preparation	Dissolve 34.66g copper sulphate tetrahydrate in water and dilute to 500ml
Usage	Reagent for reducing sugars detection
Filter paper	Filtech no: 0165, 1845

Storage and handling information

Disposal methods	16
------------------	----

Transport regulations

ERG number	154
------------	-----

