

F0876CC

Gold Line

Packaging

500ml clear plastic bottle

Physical properties, composition and data

Synonym	Alkaline potassium sodium tartrate solution; Potassium sodium tartrate alkaline solution
Liquid mixture/solution	A potassium sodium tartrate solution
Appearance	Clear, colourless, syrupy liquid
Miscibility	Miscible with water
Density (g/ml)	1.259 - 1.265
Hazardous material	Corrosive

Laboratory preparation, applications and practices

Laboratory preparation	Dissolve separately 173g potassium sodium tartrate and 50g sodium hydroxide in water and dilute when cold to 500ml Mix equal volumes of the two solutions at the time of use
Usage	Reagent for reducing sugars detection
Filter paper	Filtech no: 0225, 1839, 1893

Storage and handling information

Safety phrases	26-27-37/39-45
Risk phrases	35
Disposal methods	16

Transport regulations

Tariff code	3822.00.00
Hazardous class	8
Packing group	III
UN number	1760
ERG number	157

