

F0901NN Gold Line

powder "A"

Packaging

25g black plastic jar

F0903NN Gold Line

solution "A"

Packaging

100ml amber glass bottle

F0902NN Gold Line

powder "B"

Packaging

25g black plastic jar

F0904NN Gold Line

solution "B"

Packaging

100ml amber glass bottle

Physical properties, composition and data

Chemical group	Stains
Solid mixture	Powder: various constituents
Liquid mixture	"A" Methylene blue, Azur 1 and di-sodium hydrogen orthophosphate
Appearance	"B" Eosin yellowish, di-sodim hydrogen orthophosphate and potassium dihydrogen orthophosphate "A" Bluish purple - "B" Pinkish

Laboratory preparation, applications and practices

Usage	Detection of material parasites in thick blood smears
-------	-------------------------------------------------------

Transport regulations

Tariff code	3204.19.00
-------------	------------

