


F0910NN Platinum Line

Warranty certificate

Arsenic (As)	0.00004%
Colour Index	45350
Copper (Cu)	0.0005%
Dye content	50.5%
Iron (Fe)	0.0005%
Lead (Pb)	0.002%
Loss on drying	<10%
Water-insoluble matter	0.01%



Packaging

- 25g black plastic jar
- 100g black plastic jar
- 500g black plastic jar
- 1kg black plastic jar

Physical properties, composition and data

Chemical formula	$C_{20}H_{10}O_5Na_2$
Synonym	Fluorescein sodium; Resorcinolphthalein sodium; Sodium fluorescein; Uranin; Uranine; Uranine yellow
Chemical group	Stains
Atomic weight	376.28
Appearance	Orange red crystalline powder
Solubility	Soluble in water, slightly in alcohol
Melting point	>360°C
Boiling point	Decomposes
Density (g/ml)	Bulk: 600 kg/m ³
pH (aqueous solution)	8.3

Laboratory preparation, applications and practices

Safety phrases	25-26-39
Risk phrases	36
Usage	Analytical reagent, forensics and serology

Storage and handling information

Storage	Keep container tightly closed away from moisture and light
Disposal methods	16

Transport regulations

Tariff code	3204.12.00
-------------	------------

