


F0960HH Gold Line**Typical analysis**

Colour Index	42510
Content	96%
Loss on drying (111°C)	4.5%

Packaging

25g black plastic jar
100g black plastic jar

**Physical properties, composition and data**

Chemical formula	$C_{20}H_{20}ClN_3$
Synonym	Fuchsine; Magenta
Chemical group	Stains
Atomic weight	337.85
Appearance	Dark green powder or green crystals with bronze luster
Solubility	Soluble in water and alcohol (carmine red colour)
pH (aqueous solution)	5 - 6
Melting point	235°C with decomposition
Boiling point	Decomposes
Hazardous material	Harmful

Laboratory preparation, applications and practices

Usage	Analytical reagent. Microscopy staining
-------	---

Storage and handling information

Safety phrases	22-36/37/39-45
Risk phrases	20/21/22-40

Transport regulations

Tariff code	3204.13.00
-------------	------------

