

G0962NN Gold Line

Typical analysis

Gel strength	550-850g/ml
Loss on drying	6-16%
Particle size	78% through 28 mesh (600 micron)
	95% through 42 mesh (355 micron)
Transmittance	not less than 80%

Packaging

1 kg white plastic bucket

Physical properties, composition and data

Synonym	Mild form of Gellan gum
Appearance	White to off white powder
Boiling point	170°C
Flash point	75°C
Ignition point	490°C
pH (aqueous solution)	1% solution: 6.0 - 8.0

Laboratory preparation, applications and practices

Usage For microbiological assay media and tissue culture

Storage and handling information

Storage/handling Store in a cool dry place

Transport regulations

Tariff code 1301.90.00

