

G1005NN Gold Line

solution

Packaging

100ml amber glass bottle

500ml amber glass bottle

Physical properties, composition and data

Chemical group	Stains
Mixture	33% iodine and 67% potassium iodide
Appearance	Orange liquid
Odour	Iodine-like odour
Miscibility	Miscible with water
pH (aqueous solution)	1.005

Laboratory preparation, applications and practices

Laboratory preparation	Add 1g iodine and 2g potassium iodide to about 25ml distilled water. Make up to 300ml
Usage	Microscopy stain
Filter paper	Filtech no: 0222, 0225, 1839

Transport regulations

Tariff code	3822.00.00
-------------	------------

