


J1190NN **Gold Line**

Typical analysis

Colour Index 11050

Packaging

10g black plastic jar



Physical properties, composition and data

Chemical formula	$C_{30}H_{31}ClN_6$
Synonym	Diethylsafraninazodimethylaniline chloride
Chemical group	Stains
Atomic weight	511.09
Appearance	Black powder

Laboratory preparation, applications and practices

Usage Microscopical staining, biological stain in electrodeposition of copper and an indicator

Storage and handling information

Storage Store separated from light
 Disposal methods 3



Janus Green B indicator changes colour according to the amount of oxygen present.

Oxygen present - Indicator oxidised to blue colour

Oxygen absent - Indicator reduced to pink colour