

**K1200NN Gold Line**

heavy pure


**Typical analysis**

Complies to BP tests

Arsenic (As)	0.0002%
Calcium (Ca)	0.025%
Chloride (Cl)	0.03%
Lead (Pb)	0.002%
Loss on drying (105°C)	1.5%
Loss on ignition	15%
Sulphate (SO <sub>4</sub> )	0.1%

**Packaging**

500g white plastic jar

**Physical properties, composition and data**

Chemical formula	approx. H <sub>2</sub> Al <sub>2</sub> Si <sub>2</sub> O <sub>8</sub> ·H <sub>2</sub> O
Synonym	Bolus alba; china clay; porcelain clay; white bole; argilla
Chemical group	Silicates
Occurrence in nature	In the Earth's crust in kaolinite
Solid mixture	Hydrated aluminium silicate
Atomic weight	about 258g/mol
Appearance	White or yellowish-white, earthy mass or white powder
Density (g/ml)	1.8 - 2,6

**Laboratory preparation, applications and practices**

Usage	Analytical reagent
-------	--------------------

**Transport regulations**

Tariff code	2507.00.00
-------------	------------

