


L1270NN Platinum Line**Warranty certificate**

Assay	min 99%
Calcium (Ca)	0.0001%
Copper (Cu)	0.001%
Iron (Fe)	0.0001%
Lead (Pb)	0.0001%
Potassium (K)	0.005%

**Packaging**

100g white plastic jar
500g white plastic jar

Physical properties, composition and data

Chemical formula	LaCl ₃ ·7H ₂ O
Synonym	Lanthanum trichloride
Chemical group	Chlorides
Atomic weight	371.37
Appearance	White or transparent crystals
Solubility	Miscible with water and alcohol
Melting point	91°C
Boiling point	>1000°C
Density (g/ml)	3.842
pH (aqueous solution)	5

Laboratory preparation, applications and practices

Laboratory preparation	Treatment of lanthanum carbonate or oxides with hydrochloric acid in an atmosphere of dry hydrogen chloride
Usage	Analytical reagent

Transport regulations

Tariff code	2846.90.00
-------------	------------

