


L1290HH Gold Line**Typical analysis**

Assay	98%
Alkalis (Na + K)	1%
Iron (Fe)	0.005%

Packaging

500g white plastic jar

**Physical properties, composition and data**

Chemical formula	PbCO ₃
Synonym	Ceruse and Cerussa; Flake lead; Lead subcarbonate; White lead
Chemical group	Minerals
Occurrence in nature	In the Earth's crust in the mineral cerussite
Atomic weight	267.21
Appearance	White heavy amorphous powder
Solubility	Soluble in acids
Hazardous material	Harmful

Laboratory preparation, applications and practices

Laboratory preparation	By treating very finely divided lead in revolving wooden cylinders with dilute acetic acid and carbon dioxide
Usage	Analytical reagent

Storage and handling information

Safety phrases	53-45
Risk phrases	20/25
Disposal methods	3

Transport regulations

Tariff code	2836.70.00
-------------	------------

