

# light green sf yellowish

## L1313NN


## Gold Line

### Typical analysis

Assay	95%
C.I.	42095
Colour change	Orange to yellow
Maximum absorption	630 - 632nm

### Packaging

25g black plastic jar



### Physical properties, composition and data

Chemical formula	$C_{37}H_{34}N_2O_9S_3Na_2$
Synonyms	Lissamine green SF
Chemical group	Stains
Atomic weight	792.80
Appearance	Reddish - brown powder
Solubility	Soluble in water

### Laboratory preparation, applications and practices

Usage	Dye, biological stain
-------	-----------------------

### Storage and handling information

Storage/handling	Keep container well closed and separated from moisture and light
------------------	--

### Transport regulations

Tariff code	3204.12.00
-------------	------------

