


# lime water

L1311NN

Gold Line



## Packaging

2.5l clear plastic bottle

## Physical properties, composition and data

Chemical formula	$\text{Ca}(\text{OH})_2$
Mixture	Saturated solution of calcium hydroxide and water
Weight per litre	1.005kg
Appearance	White murky liquid
Miscibility	Miscible with water
Boiling point	100°C
Hazardous material	Corrosive

## Laboratory preparation, applications and practices

Laboratory preparation	Part calcium hydroxide and deionized water
Usage	Analytical reagent
Filter paper	Filtech no: 1285, 1803

## Storage and handling information

Storage/handling	Keep container well closed
Disposal methods	2

250