


methylated spirits

M1530FT	Silver Line
non-coloured	
Typical analysis	
Ethanol content	94%
Methanol content	6%

M1535FT	Silver Line
coloured	
Typical analysis	
Benzene (C ₆ H ₆)	1.5% v/v
Butanol	3.5% v/v
Ethanol content	94% v/v
Mineral naphtha	0.375%
Water (H ₂ O)	2.6%
Methyl violet	0.0015% w/v



Packaging
2.5l clear plastic bottle
5l clear plastic bottle
25l black plastic drum

Packaging
2.5l clear plastic bottle
5l clear plastic bottle
25l black plastic drum

Physical properties, composition and data

Chemical group	Alcohols
Synonym	Meths; Denatured alcohol
Mixture	Mixture of ethanol and methanol
Appearance	Non-coloured: clear liquid. Coloured: purple liquid
Hazardous material	Flammable and toxic

Laboratory preparation, applications and practices

Usage	Solvent. Cleaning
Filter paper	Filtech no: 8613, 8614.

Storage and handling information

Storage	Store separated from oxidants
Safety phrases	2-7-16-24
Risk phrases	11-23/25
Disposal methods	14

Transport regulations

Tariff code	2207.20.00
Hazardous class	3
Packing group	II
UN number	1230
ERG number	131

