

**N1630TC** **Gold Line**

**Packaging**

500ml amber glass bottle

**Physical properties, composition and data**

Mixture	Potassium tetraiodomercurate(II) and potassium hydroxide
Weight per litre	1.124kg
Appearance	Clear solution with yellow tinge
Density (g/ml)	1.15
pH (aqueous solution)	Strongly alkaline
Hazardous material	Toxic and corrosive

**Laboratory preparation, applications and practices**

Laboratory preparation	Solution of mercuric iodide in potassium iodide
Usage	Ammonia reaction
Filter paper	Filtech no: 0222, 0225, 1839

**Storage and handling information**

Safety phrases	26-45-36/37-61
Risk phrases	25-35-48/21/22-52/53
Disposal methods	22

**Transport regulations**

Tariff code	3822.00.00
Hazardous class	6
Packing group	III
UN number	2024
ERG number	154

