

**O1758NN Gold Line**

**Typical analysis**

Maximum absorption  
(in NaOH 0.01 mol/l) 576nm



**Packaging**

25g black plastic jar

**Physical properties, composition and data**

Chemical formula	C <sub>28</sub> H <sub>24</sub> N <sub>2</sub> O <sub>7</sub>
Synonym	Natural red 28
Chemical group	Stains/pH indicators
Mixture	Mixture of amino and hydroxyorceine and aminoorceimines
Appearance	Brownish-red microcrystalline powder
Solubility	Soluble in alcohol, acetone or acetic acid

**Laboratory preparation, applications and practices**

Laboratory preparation	By oxidation of orcinol with water in the presence of ammonia water
Usage	In microscopy for chromosome staining

**Storage and handling information**

Disposal methods	16
------------------	----

**Transport regulations**

Tariff code	3203.00.19
-------------	------------

