

P1800NN Gold Line

liquid

Typical analysis

Kinematic viscosity at 37.8°C >65cS

Packaging

500ml clear plastic bottle

2.5l clear plastic bottle

cas 8012-95-1

P1810NN Gold Line

wax refined

Typical analysis

Ash 0.05%

Packaging

500g white plastic jar

cas 80022-74-2

P1815NN Gold Line

wax pastillated

Packaging

500g white plastic jar

5Kg white plastic bucket

25Kg plastic bag

cas 80022-74-2

Physical properties, composition and data

Chemical group	Hydrocarbons
Synonym	Liquid: Petrolatum, mineral oil, paraffin oil. Solid: Hard paraffin, paraffin scale
Occurrence in nature	Wood tar
Mixture/	Mixture of hydrocarbons
Mixture	Mixture of solid hydrocarbons
Weight per litre	0.88kg
Appearance	Colourless liquid, whitish wax
Solubility/miscibility	Benzene, chloroform, ether, carbon disulphide and oils
Melting point	57 - 60°C
Density (g/ml)	0.875 - 0.905
Boiling point	350°C
Flash point	210°C
Other information	Burns with a luminous flame

Laboratory preparation, applications and practices

Usage	Analytical reagent
Filter paper	Filtech no: 0222, 0225, 1839

Storage and handling information

Disposal methods	14
------------------	----

Transport regulations

Tariff code	Liquid: 2710.00.00 Solid: 2712.20.00
-------------	---

