

# pepsin

## **P1826FF** **Gold Line**

### Typical analysis

Digestive powder 2500 units/g

### Packaging

10g white plastic jar  
100g white plastic jar  
500g white plastic jar

### Physical properties, composition and data

Synonym	Pepsinogen
Mixture	Enzyme of gastric juice
Atomic weight	34,500
Appearance	White yellowish scales or granules
Odour	Beefy smell
Solubility	Soluble in water
Incompatible substances	Alkaline substances, salts of heavy metals and tannin

### Laboratory preparation, applications and practices

Usage Controls the degradation of proteins to proteoses and peptones

### Storage and handling information

Storage Store below 4°C and separated from incompatible substances  
Disposal methods 16

### Transport regulations

Tariff code 3507.90.00

