

peptone

P1835NN

Gold Line

Packaging

100g white plastic jar

Physical properties, composition and data

Mixture	Extract from meat
Appearance	Whitish-cream powder
Odour	Beefy smell
Solubility	Soluble in water (600g/l 20°C)
pH (aqueous solution)	6 - 7

Laboratory preparation, applications and practices

Usage	Analytical reagent. Bacteriology
-------	----------------------------------

Storage and handling information

Storage	Store below 4°C
Disposal methods	16

Transport regulations

Tariff code	3504.00.00
-------------	------------

