

# potassium

## P1987CC Platinum Line

1000ppm solution

### Warranty certificate

Metal content guaranteed within 0.5% of the nominal value. Solution containing 1000ppm w/v potassium in 1N HNO<sub>3</sub>

### Packaging

500ml clear plastic bottle

### Physical properties, composition and data

Mixture	Potassium and nitric acid
Weight per litre	1.002kg
Appearance	Clear colourless liquid
Miscibility	Miscible with water and alcohol
Hazardous material	Corrosive

### Laboratory preparation, applications and practices

Usage	For spectrophotometry
Filter paper	Filtech no: 0222, 0225, 1839

### Storage and handling information

Storage	Keep container tightly closed
Safety phrases	24/25
Risk phrases	22
Disposal methods	3

### Transport regulations

Tariff code	3104.20.90
ERG number	157

