

S2270NN Gold Line

powder

Typical analysis

C.I.	50240
Dye content	>90%
Loss on drying (105°C)	<5%

Packaging

10g black plastic jar
 25g black plastic jar
 200g black plastic jar

S2271NN Gold Line

solution

Packaging

100ml amber glass bottle

**Physical properties, composition and data**

Chemical formula	$C_{20}H_{19}ClN_4$
Chemical group	Redox indicators
Atomic weight	350.86
Weight per litre	1.00kg
Appearance	Dark red-brownish powder, lustrous crystalline powder. Dark red solution
Solubility/miscibility	Water
pH (aqueous solution)	10

Laboratory preparation, applications and practices

Usage	Microscopy. Aqueous solution for bacterial spores and hegri bodies
Filter paper	Filtech no: 0222, 0225, 1839.

Storage and handling information

Safety phrases	36/37
Risk phrases	20/21/22

Transport regulations

Tariff code	Solution 3822.00.00 and powder 3204.13.00
-------------	---

