

S2290HH Gold Line**Typical analysis**

Loss on drying	2%
Water-insoluble matter	0.1%

Packaging

100g white plastic jar

Physical properties, composition and data

Chemical group	Glycosides
Synonym	Sapogenin glycoside
Mixture	Mixture of sapogenins
Appearance	Creamy yellow powder
Density (g/ml)	1.015 - 1.02
pH (aqueous solution)	4.5 - 5.5
Other information	Foams strongly when shaken with water
Hazardous material	Harmful and environment hazard

Laboratory preparation, applications and practices

Usage Analytical reagent

Storage and handling information

Safety phrases	26-37
Risk phrases	36/37
Disposal methods	16

Transport regulations

Tariff code 2938.90.00

