

schiff's reagent

S2300NN

Gold Line

Packaging

500ml amber glass bottle

Physical properties, composition and data

Chemical group	Stains
Synonym	Fuchsin-sulphurous acid
Mixture	Fuchsin, sodium bisulphite, hydrochloric acid and water
Weight per litre	1.01kg
Appearance	Purple liquid

Laboratory preparation, applications and practices

Usage	For detection of aldehydes
Filter paper	Filtech no: 0222, 0225, 1839.

Storage and handling information

Storage	Keep separated from light
---------	---------------------------

Transport regulations

Tariff code	3822.00.00
-------------	------------

