

B0241NN**Silver Line**

boiling stones

Packaging

500g white plastic jar

Physical properties, composition and data

Chemical formula	CSi or SiC
Synonym	Boiling stones
Chemical group	Carbides
Atomic weight	40.07
Appearance	Exceedingly hard, green to bluish-black, iridescent stones
Other information	2mm diameter

Laboratory preparation, applications and practices

Laboratory preparation	Heating coke and sand and salt as flux in an electric furnace
Usage	Agitation in boiling

Transport regulations

Tariff code	7018.20.00
-------------	------------

