


**S2370NN Gold Line****Typical analysis**

Assay (Argentimetric) 99.5%

**Packaging**

10g white plastic jar

**Physical properties, composition, and dataLaboratory**

Chemical formula	$\text{Ag}_2\text{CrO}_4$
Chemical group	Chromates
Atomic weight	331.77
Appearance	Dark brownish-red crystalline powder
Solubility	Soluble in nitric acid and ammonia
Density (g/ml)	5.625

**Laboratory preparation, applications and practices**

Usage Catalyst for formation of aldol from alcohol

**Transport regulations**

Tariff code 2841.50.00

