

# silver diethyldithiocarbamate

cas 1470-61-7

**S2375NN**


**Platinum Line**

## Warranty certificate

Assay (Argentimetric) min 98%

## Packaging

10g black plastic jar



## Physical properties, composition and data

Chemical formula	$(C_2H_5)_2N.CSSAg$ or $C_5H_{10}AgNS_2$
Chemical group	Carbamates
Atomic weight	256.14
Appearance	White powder
Other information	Deteriorates by natural causes

## Laboratory preparation, applications and practices

Usage	Analytical reagent
-------	--------------------

## Storage and handling information

Storage	Keep container tightly closed
Risk phrases	20/21/22
Disposal methods	13

## Transport regulations

Tariff code	2843.29.00
ECB number	216-003-7

