

S2417NN Gold Line

4-8 mesh

Typical analysis

Alkalinity (NaOH)	3.5%
Water (H ₂ O)	12 - 19%

Packaging

500g white plastic jar

Physical properties, composition and data

Mixture	Sodium hydroxide, calcium hydroxide and potassium hydroxide
Appearance	White or grey-white granules
Other information	Particle size: 4 - 8 mesh
Hazardous material	Corrosive

Laboratory preparation, applications and practices

Usage	Analytical reagent. Absorption of carbon dioxide. Drying agent
-------	--

Storage and handling information

Storage	Keep container tightly closed
Safety phrases	26-27-36/37/39-45
Risk phrases	35
Disposal methods	2

Transport regulations

Tariff code	3822.00.00
Hazardous class	8
Packing group	III
UN number	1907
ERG number	154

