

S2614NN Platinum Line**Warranty certificate**

| | |
|--------------------------------------------------|--------------|
| Assay | min 99% |
| Sensitivity to Ba, Ca, Pb and Sr | to pass test |
| Solubility in water | 0.2% |
| Suitability as indicator for sulphate titrations | to pass test |

**Packaging**

5g black plastic jar

Physical properties, composition and data

| | |
|-------------------|---------------------------------------------------------------------|
| Chemical formula | $C_6Na_2O_6$ |
| Atomic weight | 214.05 |
| Appearance | Violet crystals |
| Solubility | Soluble in water |
| Other information | Solutions are unstable and should be freshly prepared when required |

Laboratory preparation, applications and practices

| | |
|-------|---------------------------------------|
| Usage | As a reagent for barium and strontium |
|-------|---------------------------------------|

Storage and handling information

| | |
|------------------|----------------------------|
| Storage | Store separated from light |
| Disposal methods | 13 |

Transport regulations

| | |
|-------------|------------|
| Tariff code | 2942.00.00 |
|-------------|------------|