


S2716NN Gold Line**Typical analysis**

C.I.	26100
Loss on drying	<2%
Residue on ignition	<3%

Packaging

25g black plastic jar

**Physical properties, composition and data**

Chemical formula	$C_{22}H_{16}N_4O$
Synonym	Tetrazobenzene-b-naphthol; oil red; oil scarlet; solvent red 23
Chemical group	Stains
Atomic weight	352.40
Appearance	Reddish brown powder
Solubility	Chloroform and glacial acetic acid

Laboratory preparation, applications and practices

Laboratory preparation	Prepared from diazotised aminoazobenzene and 2-naphthol
Usage	Microscopy staining

Storage and handling information

Risk phrases	20/21/22
Disposal methods	16

Transport regulations

Tariff code	3204.19.00
-------------	------------

