

T2790NN Gold Line**Typical analysis**

Assay (Oxidimetric)	85%
C.I.	19140
Loss on drying	6%

Packaging

100g black plastic jar

Physical properties, composition and data

Chemical formula	$C_{16}H_9N_4Na_3O_9S_2$
Synonym	Acid yellow 23; Food yellow 4
Chemical group	Azo compounds
Atomic weight	534.37
Appearance	Bright orange-yellow powder
Solubility	Soluble in water

Laboratory preparation, applications and practices

Usage	Absorption indicator
-------	----------------------

Storage and handling information

Storage	Protect from light
Safety phrases	22-24/25-45
Risk phrases	42/43
Disposal methods	16

