

T2874TT

Gold Line

Packaging

500ml amber glass bottle

Physical properties, composition and data

Mixture

Iodine and potassium iodide

Appearance

Dark purplish-brown liquid

Incompatible substances

Alkaloids, starch and tannins

Hazardous material

Toxic

Laboratory preparation, applications and practices

Usage

Analytical reagent

Filter paper

Filttech no: 0222, 0225, 0192, 1839, 1893

Storage and handling information

Storage

Keep container tightly closed separated from incompatible substances

Transport regulations

Tariff number

3822.00.00

