

T2876CC **Platinum Line**

1000ppm solution

Warranty certificate

Metal content guaranteed within 0.5% of the nominal value. Solution containing 1000ppm w/v titanium in 2N H₂SO₄

Packaging

500ml clear plastic bottle

Physical properties, composition and data

Mixture	Titanium tetrachloride and sulphuric acid
Atomic weight	47.90
Weight per litre	1.089kg
Appearance	Colourless liquid
Hazardous material	Corrosive

Laboratory preparation, applications and practices

Usage	Spectrophotometry
Filter paper	Filtech no: 0222, 0225, 1839

Storage and handling information

Storage	Keep container tightly closed
Safety phrases	26-45
Risk phrases	34-37
Disposal methods	1

Transport regulations

Tariff code	3822.00.00
Hazardous class	8
Packing group	II
UN number	1789
ERG number	154

