

T2945NN **Platinum Line****Warranty certificate**

Assay	min 95%
Free acid	0.4%
Melting point	55°C

Packaging

100g white plastic jar
500g white plastic jar

Physical properties, composition and data

Chemical formula	$[\text{CH}_3(\text{CH}_2)_7]_3\text{PO}$ or $(\text{C}_8\text{H}_{17})_3\text{PO}$
Synonym	TOPO; Trioctylphosphine oxide
Atomic weight	386.64
Appearance	White clumpy and oily powder
Melting point	52 - 56°C

Laboratory preparation, applications and practices

Usage	Analytical reagent
-------	--------------------

Storage and handling information

Safety phrases	24/25
Risk phrases	36/38
Disposal methods	16

Transport regulations

Tariff code	2931.00.00
-------------	------------

