

**T2955NN** **Gold Line**

**Typical analysis**

Contains about 85% lactose as diluent

**Packaging**

25g white plastic jar

**Physical properties, composition and data**

Synonym	Parenzyme; Parenzymol; Tryptar; Trypure
Solid mixture	From beef pancreas
Atomic weight	24
Appearance	Yellow to grayish-yellow powder or crystals
Solubility	Soluble in water and Sorensen's sodium phosphate buffer solution
Other information	Deteriorates by natural causes

**Laboratory preparation, applications and practices**

Laboratory preparation	Proteolytic enzyme formed in the small intestine by the action of a peptidase, enterokinase, on the pancreatic cell product, trypsinogen
Usage	Enzyme

**Transport regulations**

Tariff code	3507.90.00
-------------	------------

