

<b>W2977FF</b>	<b>Silver Line</b>
low aromatic	
<b>Typical analysis</b>	
Total aromatics	4%

<b>Packaging</b>	
2.5 litre amber glass bottle	
25l metal drum	

<b>Physical properties, composition and data</b>	
Chemical group	Hydrocarbons
Mixture	A mixture of solvents, including aromatics
Weight per litre	0.78kg
Appearance	Clear, colourless liquid
Flash point	39°C
Density (g/ml)	0.75 - 0.78
Hazardous material	Flammable

<b>Laboratory preparation, applications and practices</b>	
Usage	Analytical reagent
Filter paper	Filtech no: 0222, 0225, 1839

<b>Storage and handling information</b>	
Storage	Store in a fireproof place
Safety phrases	24/25
Risk phrases	10-20/21/22

<b>Transport regulations</b>	
Hazardous class	3
Packing group	II
UN number	1300
ERG number	128

