

W2980NN**Gold Line**

powder

Packaging

25g black plastic jar

W2985FT**Gold Line**

solution

Packaging

100ml amber glass bottle

500ml black plastic bottle

Physical properties, composition and data

Chemical group	Stains
Mixture	Eosin-methylene blue
Appearance	Greenish-blue liquid. Black powder
Odour	Musty odour
Flash point	12°C
Hazardous material	Solution: flammable and toxic

Laboratory preparation, applications and practices

Usage	Microscopy
Filter paper	Filtech no: 0222, 0225, 1839

Storage and handling information

Safety phrases	7-16-45
Risk phrases	11-23/25
Disposal	15

Transport regulations

Tariff code	3212.90.00
Hazardous class	Solution: 3
Packing group	Solution: II
UN number	1992
ERG number	Solution: 131

